Healthy water, healthy life.





Filter Performance

Verified by independent laboratories

Candle Grade 49mmø Selection of the most appropriate grade of candle allows filter performance to be matched to the requirements of the local water conditions.		Sterasyl® Silver impregnated ceramic microfilter	Super Sterasyl® Sterasyl® microfilter plus granular activated carbon	ATC Super Sterasyl® Sterasyl® microfilter plus granular activated carbon and lead removal media	Supercarb® Sterasyl® microfilter plus carbon block core	Chlorasyl® Sterasyl® microfilter plus carbon block core	Ultracarb® Sterasyl® microfilter plus carbon block core and lead removal media
	Absolute (defined as >99.99%) Nominal (defined as >99.9%)	0.9 micron 0.5 to 0.8 micron	0.9 micron 0.5 to 0.8 micron	0.9 micron 0.5 to 0.8 micron	0.9 micron 0.5 to 0.8 micron	0.9 micron 0.5 to 0.8 micron	0.9 micron 0.5 to 0.8 micron
Working Pressure (for pressure filte	er use)						
	Minimum Maximum	10 psi 125 psi	10 psi 125 psi	N/A N/A	10 psi 125 psi	10 psi 125 psi	10 psi 125 psi
Working Parameters							
	Working temperature range Working ph range Suitable for use in Gravity filters Recommended change frequency	5 - 30 (°C) 5.5 - 9.5 Yes 12 months	5 - 30 (°C) 5.5 - 9.5 Yes 6 months	5 - 30 (°C) 5.5 - 9.5 Yes 6 months	5 - 30 (°C) 5.5 - 9.5 No 6 months	5 - 30 (°C) 5.5 - 9.5 No 6 months	5 - 30 (°C) 5.5 - 9.5 No 6 months
Flow Rate							
Unrestricted Flow at 3 Bar Pressure To achieve maximum performance	Litres per minute US gallons per minute Litres per minute	5 1.33 N/A	4.5 1.2 1.5	N/A N/A 1.2 I/hr under gravity	3.7 1 1.9	3.7 1 1.9	3.3 0.9 1.9
	US gallons per minute	N/A	0.4	0.3 g/hr under gravity	0.5	0.5	0.5
Capacity Before replacement to quarantee performance	Litres US gallons	10,000 2600	2000 535	1500 400	3800 1000	3800 1000	2300 600
Quality Approval	S certified version available WRAS approved Turbidity reduction to NSF std. 53	YES+ YES >98%	NO YES >98%	NO NO >98%	YES+ YES >98%	No No >98%	YES+ YES >98%
Pathogenic Organisms							
% Bacteria Removal	E. Coli / Cholera / Shigella / Typhoid / Klebsiella Terrigena	>99.99%	>99.99%	>99.99%	>99.99%	>99.99%	>99.99%
% Cyst Removal	Cryptosporidium Giardia	>99.99% >99.99%	>99.99% >99.99%	>99.99% >99.99%	>99.99% >99.99%	>99.99% >99.99%	>99.99% >99.99%
Trace Organics Removal							
Insecticides Herbicides Phenols Polyaromatic Hydrocarbons Trihalomethanes	Lindane @ 0.1ppb presence Atrazine @ 1.2ppb presence TCP @ 1.2ppb presence PAHs @ 0.2ppb presence Chloroform @ 150ppb presence	N/A N/A N/A N/A	NO DATA NO DATA NO DATA NO DATA NO DATA	NO DATA NO DATA NO DATA NO DATA NO DATA	>85% >85% >50% >95% >50%	NO DATA NO DATA NO DATA NO DATA NO DATA	>85% >85% >50% >95% >50%
	онногоготт — гоорри ргезепсе	IWA.	NO DAIA	NO DAIA	~00 /0	NO DAIA	Z0070
Inorganics Removal Free Chlorine Removal	2mg cl/l challenge	N/A	Under Gravity >95% Under Pressure >50%	>95% under gravity	>97% NSF Class 1 std. 42	>97% (1000 USG)	>96% NSF Class 1 std. 42
Monochloramine Reduction Tested to 500 US gallons (1893 litres Tested to 700 US gallons (2650 litres	, 01	N/A N/A	N/A N/A	N/A N/A	N/A N/A	95% 95%	N/A N/A
Lead Removal	@ pH 6.5 after 2,244L @ pH 8.5 after 2,244L	N/A N/A	N/A N/A	>94% under gravity N/A	N/A N/A	No No	>99.3% NSF std. 53

Scale Inhibition

Slow release scale inhibition available with Ultracarb SI® and Supercarb SI® ceramic filter candles.

State restrictions may limit the availability of certain products in California, Iowa and Wisconsin. For product availability in these States please refer to our website: www.faireyceramics.com



⁺ NSF certified replacement elements are Doulton Sterasyl®, Doulton® Supercarb® and Doulton® Ultracarb® as verified by testing in HIP Inline systems.